

Increasing Selectivity and Sensitivity in Ion Mobility Spectrometry using Ion Beam Modulation

M. Todd Griffin; Jack E. Fulton, Jr.; Robert F. McAtee;

Chemical, Biological, and Explosives Detection Group NSWC, Crane Division Crane, IN 47522

<u>Distribution Statement A</u> - Approved for public release; distribution unlimited.



maintaining the data needed, and c including suggestions for reducing	election of information is estimated to completing and reviewing the collect this burden, to Washington Headquuld be aware that notwithstanding an OMB control number.	ion of information. Send comments arters Services, Directorate for Infor	regarding this burden estimate of mation Operations and Reports	or any other aspect of the 1215 Jefferson Davis	nis collection of information, Highway, Suite 1204, Arlington	
1. REPORT DATE 01 OCT 2005		2. REPORT TYPE N/A		3. DATES COVE	RED	
4. TITLE AND SUBTITLE					5a. CONTRACT NUMBER	
Increasing Selectivity and Sensitivity in Ion Mobility Spectrometry using Ion Beam Modulation					5b. GRANT NUMBER	
Ion Deam Modulation				5c. PROGRAM ELEMENT NUMBER		
6. AUTHOR(S)				5d. PROJECT NUMBER		
					5e. TASK NUMBER	
					5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Chemical, Biological, and Explosives Detection Group NSWC, Crane Division Crane, IN 47522					8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)		
		11. SPONSOR/MONITOR'S REPORT NUMBER(S)				
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release, distribution unlimited						
	otes 51, Proceedings of t Research, 17-20 No					
14. ABSTRACT						
15. SUBJECT TERMS						
16. SECURITY CLASSIFIC	17. LIMITATION OF	18. NUMBER	19a. NAME OF			
a. REPORT unclassified	b. ABSTRACT unclassified	c. THIS PAGE unclassified	ABSTRACT UU	OF PAGES 18	RESPONSIBLE PERSON	

Report Documentation Page

Form Approved OMB No. 0704-0188



Chem/Bio/Explosive Detection



FACILITIES

- Unique 25,000 Sq. Ft. Chemical/Biological Detection Center
- Physical Plant Replacement Value of \$4.0M
- Specialized Chemical/Biological Detection Equipment with Replacement Value of \$5.0M
 - Thermal Collimators
 - Chemical Vapor Diagnostic Test Sets
 - Chemical Vapor Generators
 - Chemical Hazard Containment Room

CAPABILITIES

- Develop & Review Specifications
- Develop Procurement Technical Data Packages
- Contract for Products, Components or Services
- Acceptance Test New Systems
- Repair Chemical Detection Devices
- Provide On-Site Fleet Support and Training
- Initiate & Implement Engineering Change Proposals
- Provide Sustainment Logistics Support for New & Existing Systems
- Perform Failure Analysis
- Hold Navy Radioactive Material Handling Permit

PERSONNEL

- Highly Technical Workforce of 51
 Professionals with Extensive
 Chemical/Biological/Explosive
 Detection Experience
- Skills Mix of Thermal Imagers, and IMS Expertise
- 19 Engineers and Scientists
- 32 Skilled Technicians/Logisticians/ Equipment Specialists



Code 805D Customers

Navy

- Naval Sea Systems Command
- Naval Special Warfare Command
- Navy Inventory Control Point, Mechanicsburg
- Navy Explosive Ordnance Disposal Units
- Navy Construction Battalions
- Office of Naval Research

Other

- Defense Threat Reduction Agency
- Soldier, Biological, Chemical Command
- Joint Program Office Biological Defense





Code 805D Tasking

- Acquisition Engineering Support
- Production Acceptance Testing
- Installation & Fielding Support
- In-Service Engineering Support
 - Product Improvement Programs
 - Configuration Management- Hardware and Software
 - Intensified Leak Test Tracking for Radioactive Sources
- Depot Maintenance and Repair
- Support New Construction Ships
- Product Integration (System Design and Integration)





IMS in the Field



Automatic Chemical Agent Detection Alarm (ACADA)



Improved Chemical Agent
Monitor (ICAM)



Improved Point Detection Sys. (IPDS)








VaporTracer II (Explosives)






Ion Mobility Spectrometry: System Block Diagram

Atmosphere






Field Portable IMS System





Typical IMS Signature: Positive Cell




DPM from Confidence Tube





IMS Response at Low Concentration




 0.027 mg/m^3





BNG Driver Timing



50% Duty Cycle






Background

11 of 18


Wavelet Decomposition





DPM


Wavelet Decomposition








Time-of-flight Measurement







Wavelet Decomposition




14 of 18

Harnessing the Power of Technology for the Warfighter






Neural Networks



Adaptive Learning Rule:

$$W^{new} = W^{old} + (t - a) \bullet X^{T}$$



- Benefit of Multi-layered NN: They can approximate any function
- <u>Difficulty</u>: Updating weights in the hidden layers
- <u>Solution</u>: Back-propagation of error
 - Gradient Descent:

$$\Delta W = -\eta \frac{\partial E}{\partial W}$$

$$E = \frac{1}{2} \sum_{i} (t_{i} - a_{i})^{2}$$

$$a_1 = f(x_1 w_1 + x_2 w_2 + \Lambda + x_n w_n) = f\left(\sum_{i=1}^n x_i w_i\right) = f(W^t T)$$



Results

- Collected data sets using ion beam modulation
 - Background
 - DPM
 - 1-propanol
- Performed Wavelet Decomposition
- Used 5 new statistical evaluators
- Trained 5-node Neural Network (Feed forward)

Analyte	Results	% Correct
Background	3/3	100%
DPM	5/5	100%
Propanol	3/3	100%



Conclusions and Future Work

- Wavelets can help improve S/N
- Wavelets can be used to decompose signal for further analysis
- Ion beam modulation improves selectivity and sensitivity
- Neural networks can be used to detect modulated ion beam data

- Investigate modulation schemes
- Determine algorithm robustness
- Decrease IBM scan time
- Decrease IBM data acquisition rate
- Expand training sets
- Look at other metrics (cell temperature, humidity)





Acknowledgements

Funding:

Integrated Detection of Hazardous Materials (IDHM) Program, a Department of Defense project managed jointly by Center for Sensing Science and Technology, Purdue University, and Naval Surface Warfare Center, Crane, Indiana.

Collaborators:

Rong Gao*, Lefteri H. Tsoukalas Applied Intelligent Sensing Lab (AISL) School of Nuclear Engineering, Purdue University West Lafayette, IN 47907

